

250 mL 25x Wash Solution Concentrate for KIT 029

PRODUCT NO. **WASH 029-250**

PRESENTATION Preparation: 25x Wash Solution Concentrate
Content: 250 mL
Storage: Store at 2-8°C. Do not freeze.
Stability, undiluted: See label
Stability after dilution: 4 weeks at 2-8°C. Do not freeze.

PREPARATION Calculate the amount of 25x Wash Solution Concentrate needed for producing the required volume of diluted wash solution and add distilled or deionized water in an 1/25 ratio. Mix thoroughly and store at 2-8°C after use.

INTENDED USE The 250 mL 25x Wash Solution Concentrate for KIT 029 is intended for use with the MBL Oligomer ELISA Kit when used with automated ELISA washing systems. The kit already contains 30 mL 25x Wash Solution Concentrate, which should be sufficient for manual washing.

INSTRUCTIONS Please refer to the instructions for use for KIT 029 MBL Oligomer ELISA Kit for further information.

CONDITIONS

All products are supplied on the understanding that they are for in vitro use only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control