

**Anti Gc-globulin (human)
Mouse Monoclonal Antibody**Subclass: IgG₁/κ

PRODUCT NO.	HYB 249-02																
PRESENTATION	Preparation: Protein-A/G purified Content: 1 ml, 1 mg/mL Solvent: 0.01 M phosphate buffer, pH 7.4, with 0.5 M NaCl and 15mM sodium azide Storage: In the dark at 4-8°C																
ANTIGEN	Gc-globulin is a plasma protein produced in the liver. Amongst its ligands are vitamin D, thus Gc-globulin is also called vitamin D-binding protein. Gc-globulin is furthermore part of the actin scavenging system, binding and removing monomeric actin from the blood stream. The molecular mass of Gc-globulin is approximately 50 kDa. The concentration of Gc-globulin in human plasma is app. 400 µg/ml (1).																
IMMUNOGEN	Gc-globulin isolated from human plasma																
SPECIFICITY	HYB 249-02 is specific for human Gc-globulin.																
EPI TOPE SPECIFICITY	Epitope specificity differs from that of HYB 249-01, HYB 249-05 and HYB 249-10.																
REACTIVITY	HYB 249-02 reacts strongly with Gc-globulin. Strong reaction is seen in ELISA with Gc-globulin coated directly onto the microtiter well and when tested in sandwich ELISA in combination with a polyclonal antibody against Gc-globulin (e.g. DAKO A 0021). In Western blotting after SDS-PAGE, HYB 249-02 reacts with Gc-globulin in both reduced as well as unreduced forms.																
CULTURE MEDIUM	RPMI 1640 with 10% fetal calf serum																
FUSION PARTNER	X63-Ag8.653.																
IMMUNIZATION	Female CF1 x BALB/c mice, immunized i.p. with immunogen adsorbed onto Al(OH) ₃																
APPLICATION	<table border="1"> <thead> <tr> <th>Method</th> <th>Usability</th> <th>Dilution guideline</th> <th>References</th> </tr> </thead> <tbody> <tr> <td>ELISA</td> <td>Yes</td> <td>1:20,000</td> <td></td> </tr> <tr> <td>Immunoblotting</td> <td>Yes</td> <td></td> <td></td> </tr> <tr> <td>Immunohistochemistry</td> <td>Not determined</td> <td></td> <td></td> </tr> </tbody> </table> <p>The dilution guideline for ELISA is based on sandwich ELISA in combination with a polyclonal antibody against the antigen. Users should determine the optimal dilutions for their own purpose.</p>	Method	Usability	Dilution guideline	References	ELISA	Yes	1:20,000		Immunoblotting	Yes			Immunohistochemistry	Not determined		
Method	Usability	Dilution guideline	References														
ELISA	Yes	1:20,000															
Immunoblotting	Yes																
Immunohistochemistry	Not determined																
REFERENCES	1. Masuda S, Okano T, Osawa K, Shinjo M, Suematsu T, Kobayashi T (1989) Concentrations of vitamin D-binding protein and vitamin D metabolites in plasma of patients with liver cirrhosis. J Nutr Sci Vitaminol 35:225-34.																

CONDITIONS

All products are supplied on the understanding that they are for in vitro use only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The animals from which this product was derived have not been exposed to or inoculated with any livestock or poultry disease agents exotic to the United States or Western Europe, and did not originate from facilities where work with exotic disease agents affecting livestock or avian species is carried out.