

**Anti Gc-globulin (human)
mouse monoclonal antibody**

Subclass: IgG2b/k

PRODUCT NO.

HYB 249-05

PRESENTATION

Preparation: Protein-A/G purified
 Content: Available in 200 µL and 1 mL, 1 mg/mL
 Solvent: 0.01 M phosphate buffer, pH 7.4, with 0.5 M NaCl and 15 mM sodium azide
 Storage: In the dark at 4-8°C

ANTIGEN

Gc-globulin is a plasma protein produced in the liver. Amongst its ligands are vitamin D, thus Gc-globulin is also called vitamin D-binding protein. Gc-globulin is furthermore part of the actin scavenging system, binding and removing monomeric actin from the blood stream. The molecular mass of Gc-globulin is approximately 50 kDa. The concentration of Gc-globulin in human plasma is app. 400 µg/ml (1).

IMMUNOGEN

Gc-globulin isolated from human plasma.

SPECIFICITY

HYB 249-05 is specific for human Gc-globulin.

EPI TOPE SPECIFICITY

Epitope differs from that of HYB 249-01, HYB 249-02 and HYB 249-10.

REACTIVITY

HYB 249-05 reacts strongly with Gc-globulin. Strong reaction is seen when tested in sandwich ELISA in combination with a polyclonal antibody against Gc-globulin (e.g. DAKO A0021).

CULTURE MEDIUM

RPMI 1640 with 10% fetal calf serum

FUSION PARTNER

X63.Ag8.653

IMMUNIZATION

Female CF1 x Balb/c mice, immunised i.p. with antigen adsorbed onto Al(OH)₃.

APPLICATION

Method	Usability	Dilution guideline	References
ELISA	Yes	1:25,000	
Immunoblotting	Not determined		
Immunohistochemistry	Not determined		

The dilution guideline for ELISA is based on sandwich ELISA in combination with a polyclonal antibody against the antigen. Users should determine the optimal dilutions for their own purpose.

REFERENCES

1. Masuda S, Okano T, Osawa K, Shinjo M, Suematsu T, Kobayashi T (1989) Concentrations of vitamin D-binding protein and vitamin D metabolites in plasma of patients with liver cirrhosis. J Nutr Sci Vitaminol 35:225-34.

CONDITIONS

All products are supplied on the understanding that they are for in vitro use only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The animals from which this product was derived have not been exposed to or inoculated with any livestock or poultry disease agents exotic to the United States or Western Europe, and did not originate from facilities where work with exotic disease agents affecting livestock or avian species is carried out.