

**Anti-Complement component C2 (human)
 Mouse monoclonal antibody**

Subclass: IgG_{2a}/k

PRODUCT NO.

HYB 050-04

Clone: 5C3

PRESENTATION

Preparation: Protein-A/G purified

Content: Available in 200 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

Solvent: 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

Storage: 4-8°C without exposure to light. No precautions necessary during handling.

ANTIGEN

Complement component C2 is a component of the classical complement pathway and is involved in activation of C3 and C5. C2 consists of a single polypeptide chain with a molecular mass of 102 kDa (2).

IMMUNOGEN

Complement component C2 isolated from human plasma adsorbed onto aluminum hydroxide gel (1)

SPECIFICITY

HYB 050-04 is specific for C2 from human plasma/serum. Cross-reactivity to C2 from other species has not been tested.

EPI TOPE SPECIFICITY

Epitope specificity differs from that of HYB 050-05 and HYB 050-08 as determined by inhibition ELISA.

REACTIVITY

HYB 050-04 reacts with human C2 in both native and denatured forms. In Western blotting after SDS-PAGE, HYB 050-04 reacts with C2 and a subfraction of C2 believed to be C2b, both under reducing and non-reducing conditions. In ELISA, strong reaction is seen with C2 directly coated onto the microtiter well.

CULTURE MEDIUM

RPMI 1640 with 10% fetal calf serum

FUSION PARTNER

X63-Ag8.653.

IMMUNIZATION

Female CF1 x BALB/c mice immunized by intraperitoneal injection

APPLICATION

Method	Usability	References
ELISA	Yes	1
Immunoblotting	Yes	1
Immunohistochemistry	Not determined	

REFERENCES

1. Stenbaek EI, Koch C, Barkholt V, Welinder KG (1986) Human complement component C2: production and characterization of polyclonal and monoclonal antibodies against C2. Mol Immunol 23:879-886.
2. Law SKA, Reid KBM (1988) Complement. In: In Focus (Ed. Male D) IRL Press: Oxford.

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.